
	Minivisc® plus
Intrinsic Viscosity	1.5 - 4.5 m ³ /kg
pH	6.8 - 7.6
Endotoxin	NMT 0.2 EU/ml
Protein residuals	0.1%
Refractive index	1.330 - 1.350
Particles ≥10 µm	≥ 3000 per container
≥25 µm	≥ 300 per container
Shear Viscosity app.	323 Pa.s

MINIVISC® PLUS

Sodium hyaluronat 1,4%

1050
SCANDINAVIAN HA

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E-mail: info@bohusbiotech.com
Web: www.bohusbiotech.com

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

CH

REP

INSTRUCTIONS FOR USE	
Minivisc® plus	
Formula	
Each ml contains	Sodium hyaluronate 1.4%
Sodium hyaluronate	14.0 mg
Disodium phosphate dihydrate	1.4 mg
Sodium chloride	8.3 mg
Potassium dihydrogen phosphate	0.26 mg
Water for injection	q.s
Osmolality	200-400 mOsm/kg

Description
Minivisc® plus is an ophthalmic viscoelastic device (OVD) to be used as a surgical aid in ophthalmic surgery. The device is a sterile, nonpyrogenic, viscoelastic solution of sodium hyaluronate dissolved in a physiological buffer at pH 6.8-7.6. It is a clear solution supplied in a disposable glass syringe. A sterile single-use cannula is enclosed in each box.

Characteristics
Hyaluronic acid is a native substance with high biocompatibility. The most important properties of Hyaluronic acid are to protect, lubricate and support delicate cells and tissues. Minivisc® plus forms a thin layer without causing compression of delicate cells or tissues.

The hyaluronic acid for the product Minivisc® plus is produced by bio fermentation from Streptococcus bacteria.

Intended use
Used temporary during intraocular surgery for protection, lubrication and support of delicate cells or tissues in ophthalmic procedures. For intraocular use only, in Ophthalmic procedures including:

1. Cataract surgery and IOL implantation
2. Glaucoma surgery
3. Corneal transplantation
4. Anterior segment surgery

Contraindications
None known at present.

Precautions
• A transient increase in postoperative intraocular pressure (IOP) has been reported with the use of viscoelastic preparation. The increased pressure will reach a maximum four to eight hours following surgery. The IOP should be checked the days following operation.

- Residual viscoelastic not removed at the end of the surgical procedure is flushed out naturally via the trabecular system and Schlemm's canal. However, there is the risk of blockage of the drainage channels, which can lead to an increase in intraocular pressure. To prevent the risk of IOP, irrigation /aspiration is used at the end of the surgical procedure to remove residual viscoelastic. In case of IOP increase, treat with the preferred therapy. Patients who have been diagnostic with wide-angle glaucoma, severe myopia, diabetic retinopathy or uveitis are at particular risk of increased intraocular pressure.
- There is no evidence concerning the safety of Minivisc® plus in human pregnancy and lactation. Administration during pregnancy and lactation is at the discretion of the ophthalmic surgeon.

Interactions
The user should be aware of the possibility that drugs applied or injected at the same time and same place as Minivisc® plus may have a longer residue time due to the matrix structure of Minivisc® plus. Minivisc® plus becomes cloudy if mixed with benzalkonium chloride, chlorhexidine or certain washing detergents. If cloudiness occurs, the solution must be removed from the eye.

The product is intended for intraocular use only, as an aid in eye surgery. It is instilled intraocularly by means of the enclosed thin, bended cannula. The products may only be used by medical professionals with adequate training.

Minivisc® is produced from bacterial Sodium Hyaluronate. Since it contains very small amounts of potentially immunogenic material, a possibility exists that local reactions may appear. To minimize such a reaction the patients should be asked whether they are allergic to bacterial Sodium Hyaluronate.

- Warnings**
- At the end of the operation the solution shall be removed by irrigation or aspiration.
 - The syringe, cannula and any unused material must be discarded after the treatment. Disposal should be in accordance with accepted medical practice and applicable national or local guidelines.
 - The product may only use by medical professionals with adequate training.
 - This product is for single use only. If used for repeated administrations, sterility can not be guaranteed.
 - Do not reserialize the prefilled syringe and cannula.
 - Do not use if package is damaged or opened.
 - Do not use after expire date printed on the product.

How Supplied
Minivisc® plus is a clear solution supplied in a disposable glass syringe. A sterile single-use 27G cannula is enclosed in each box. Each syringe contains 1.0 ml of 1.4% of sodium hyaluronate.

- Instruction for use**
- Before use allow the preparation to equilibrate to room temperature for 20 minutes.
 - Remove the tip cap.
 - Mount the cannula by using the thumb and forefinger to hold firmly around both the glass syringe barrel and the luer-lock adapter.
 - Grasp the lower part of needle (luer lock) with the other hand. To facilitate proper assembly, both push and rotate firmly.

Shelf life is 3 years from manufacturing date.

SYMBOLS	
	Consult instructions for use
	Do not use if package material is damaged
	Keep away from sunlight
	Single-use only
	Steam sterilized
	Sterilized using ethylene oxide
	Do not resterilize
	Store at 2-25° C. Protect from light and freezing
	Keep dry
	Manufacturer
	Expiry date (year/month)
	Date of manufacture (year/month)
	Version number
	Do not use if package material is damaged
	Keep away from sunlight
	Single-use only
	Steam sterilized
	Sterilized using ethylene oxide
	Do not resterilize
	Store at 2-25° C. Protect from light and freezing
	Keep dry
	Manufacturer
	Expiry date (year/month)
	Date of manufacture (year/month)
	Version number

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