
DERMAFILL® GLOBAL XTRA VOLUME ULTRA S-INTENSE LIPS

VER 2024-06 1053

DERMAFILL is a SBS-MED brand

Distributed by
Infinity Cosmetics SAS
ZA Les Grands Champs, 44300 Nantes - France
www.dermafill.eu

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DESCRIPTION

DERMAFILL® are transparent, injectable implants in the form of a viscous, sterile, non-pyrogenic gel. **DERMAFILL®** are composed of cross-linked sodium hyaluronate of non-animal origin.

DERMAFILL® are supplied in a graded, disposable, 1.0 ml glass syringe. The product is for single use only. Two disposable sterile needles are provided with each syringe. Information about the needle size and sterilization method can be found on the outer package of the needle. The volume in each syringe is stated on the outer package. Two patient record labels will be included in each package to ensure product traceability.

INDICATIONS

DERMAFILL® Global Xtra is intended for treatment of fine lines, superficial wrinkles and lip correction. The product should be injected intradermally.

DERMAFILL® Lips is intended to restore and/or augment lip volume. The product should be injected into the submucosal layer.

DERMAFILL® Volume Ultra and **DERMAFILL® S-Intense** are intended for treatment of medium to deep wrinkles and folds, and lip enhancement. The product is to be injected intradermally.

WARNINGS AND PRECAUTIONS

- Always observe the general precautions applicable to intradermal injections.
- **DERMAFILL®** contain traces of BDDE (1,4-butanediol diglycidyl ether, less than 0.0001%) used as cross-linking agent, which in rare cases can cause irritation to the dermis.
- **DERMAFILL®** may only be used as an intradermal implant by properly trained, qualified physicians or the authorized staff members in accordance with national legislation.
- The customary intradermal injection technique must be used.
- Ensure that the integrity of the syringe and needle has not been compromised and check the expiry date prior to use.
- **DERMAFILL®** with other products or inject the products into the blood stream.
- Do not inject **DERMAFILL®** directly adjacent to a permanent implant or in connection with other on-going or planned skin treatments, e.g. laser treatment, in the same area.
- Do not use **DERMAFILL®** on patients with infections or inflammations in the intended injection area.
- Please note that injections of **DERMAFILL®** can reactivate a latent herpes simplex virus or cause angioedema-type swelling.
- If using substances that may have an effect on the platelet function, such as anti-inflammatory drugs or aspirin, the patient may to a larger extent experience bruising or bleeding at injection sites.
- **DERMAFILL®** have not been tested on pregnant or breast-feeding women or on children.

CONTRAINDICATIONS

- **DERMAFILL®** must not be administered to patients with any acute or chronic skin disease or inflammation (such as pimples, rashes or hives) within or close to the area selected for correction.
- **DERMAFILL®** must not be administered to patients with bleeding disorders or in patients who are taking thrombolytics or anticoagulants.

PATIENT INFORMATION

- After treatment with **DERMAFILL®**, undesirable short-term reactions may occur at the injection site, such as pain, tenderness, redness or swelling. The reactions normally appear shortly after treatment and may last up to 7 days for injections into the skin and up to 14 days for injection into the lips. If the patient experiences severe discomfort or if the site gives rise to other conditions, such as oedema or prolonged swelling, contact the physician who administered the treatment.
- **DERMAFILL®** break down naturally over time and are reabsorbed by the surrounding tissue. Durability varies with the patient and the site of the implant.
- The patient must not expose the treated area to extreme temperatures (e.g. cold, heat, sun bathing, solarium etc.) during the first week after treatment or until any undesirable reactions have fully healed. Strenuous exercise should be avoided within the first day or two following application.
- The patient should be informed not to apply make-up prior to, and up to at least 12 hours after, treatment.

DIRECTIONS FOR USE

N.B. DO NOT USE THE PRODUCT IF THE PACKAGE HAS BEEN OPENED OR DAMAGED!

General directions:

- Strict aseptic technique must be performed at all times.
- The products are to be injected intradermally into the dermis under the wrinkle or fold, increasing volume and adding lift and shape. When correcting wrinkles and folds, although they will be smoothed out, the line will always remain.
- The result of the treatment is based on using a correct injection technique **DERMAFILL®** must therefore only be administered by authorized personnel in accordance with local legislation.
- The expected results of treatment must be assessed for each patient individually. Patients with unrealizable expectations should not be subject to treatment.
- The patient's medical background should always be checked before treatment in order to ensure that there are no contraindications. The patient should be informed about the indications, contraindications as well as warnings and precautions.
- The need for any pain relief should be determined.
- In order to achieve and sustain the maximum desired result, the patient should be informed of the need for additional treatment(s).

Directions before and during treatment:

- For information on how to assemble the needle and syringe, please see below.
- To avoid the risk of infection during implant, ensure that the skin has been completely disinfected with alcohol, or other antiseptic solution. Standard practice to prevent infections is always to be exercised.
- Do not inject the products too deep as this will shorten the effect of the application.
- Do not inject the products too superficially as this may result in visible lumps and/or discoloration of the skin.
- Inject slowly into the dermis while pulling the needle backwards.
- Complete the injection a couple of millimetres before the entry point in the skin to prevent material from leaking out of the injection site.
- Correction should only be made to 100% of the desired volume effect. Do not overcorrect.
- It is recommended to change needle/cannula for each new treatment site.
- For each treatment site a maximum dosage of 2 ml per treatment session is recommended.
- Following an injection of **DERMAFILL®**, the skin may be lightly massaged to encourage the product to follow the contours of the surrounding tissue.
- If stored below room temperature, allow product to equilibrate to room temperature for at least 20 minutes prior to use.

Disposal and single use:

- The syringe, needle and any unused material must be discarded immediately after the treatment.
- No material remaining after treatment may be used for more than one patient, be saved or be re-sterilized.
- Disposal should be in accordance with accepted medical practice and applicable national or local laws and guidelines.

N.B. ANY ADVERSE EVENTS MUST BE REPORTED TO THE LOCAL DISTRIBUTOR OR DIRECTLY TO THE MANUFACTURER.

ASSEMBLY OF NEEDLE AND SYRINGE

To avoid separation of the needle and syringe during implant, please notice the following instructions:

- Remove the tip cap as displayed (pictures 1 and 2).
- Hold the syringe by using the thumb and forefinger (3) to hold firmly around both the glass syringe barrel and the luer-lock adapter, then grasp the needle, (still inside the needle protector) with the other hand (3).
- To facilitate proper assembly, both push and rotate firmly at the same time (4).
- Remove the needle protector in a straight, forward motion (5).

Strict aseptic technique must be performed at all times.

STORAGE AND SHELF LIFE

- Store the products at a temperature of 2-25°C (36-77°F).
- The products must not be frozen.
- Keep out of direct sunlight.
- The expiry date of the product is indicated on the package.

SYMBOLS

- Consult instructions for use
- Single-use only
- Do not resterilize
- Steam sterilized product
- Sterile / EO needle
- Store at 2-25° C
- Keep dry
- Keep away from sunlight
- Do not use if package is damaged
- Non pyrogenic
- Caution
- Medical device name
- Cannula
- Syringe
- Version number
- Batch code
- Use by date
- Manufacturer
- CE-mark of conformity (0483)

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