30-29396
VER 2021-03

0483
0197

Manufactured by:

Bohus BioTech AB
Trädgårdsgatan 4
SE-452 31 Strömstad, Sweden
E-mail: info@bohusbiotech.com
Web: www.bohusbiotech.com

Distributed by:

Qatuainah Medical Group
KIPCO Tower, Floor 45, Office 2,
Al Shuhada St., Kuwait City
Email: info@qmg-kuwait.com
Web: www.qmg-kuwait.com

Elevate® KISS
Cross-linked sodium hyaluronate 20 mg/ml

Elevate® ULTRA
Cross-linked sodium hyaluronate 20 mg/ml

Elevate® VOLUMA
Cross-linked sodium hyaluronate 20 mg/ml

en INSTRUCTIONS FOR USE

Elevate® Kiss/Elevate® Ultra/Elevate® Voluma

PREPARATION

Contents	Elevate Kiss	Elevate Ultra	Elevate Voluma
Sodium hyaluronate, cross-linked	20 mg/ml	20 mg/ml	20 mg/ml
Sodium chloride solution	q.s	q.s	q.s
Physiological buffer	q.s	q.s	q.s

DESCRIPTION

Elevate Kiss, Elevate Ultra and Elevate Voluma are transparent, injectable implants in the form of a viscous, sterile, non-pyrogenic gel.
Elevate Kiss, Elevate Ultra and Elevate Voluma are composed of cross-linked sodium hyaluronate of non-animal origin.
Elevate Kiss, Elevate Ultra and Elevate Voluma are supplied in a graded, disposable, 1.0 ml glass syringe. The products are for single use only. Two disposable sterile needles are provided with each syringe. Information about the needle size and sterilization method can be found on the outer package of the needle. The volume in each syringe is stated on the outer package. Patient record labels are included in each package to ensure product traceability.

INDICATIONS

Elevate Kiss is intended for treatment of fine lines, superficial wrinkles and lip correction. The product should be injected intradermally.
Elevate Ultra is intended for treatment of medium to deep wrinkles and folds, and lip enhancement. The product should be injected intradermally.
Elevate Voluma is intended for facial augmentation, including cheek or chin augmentation, to correct age-related volume deficit.

WARNINGS AND PRECAUTIONS

- Always observe the general precautions applicable to intradermal injections.
- Elevate Kiss, Elevate Ultra and Elevate Voluma contain traces of BDDE (1,4-butanediol diglycidyl ether, less than 0.0001%) used as cross-linking agent, which in rare cases can cause irritation to the dermis.
- Elevate Kiss, Elevate Ultra and Elevate Voluma may only be used as an intradermal implant by properly trained, qualified physicians or aesthetic staff members in accordance with local regulations.
- The customary intradermal injection technique must be used during treatment.
- Ensure that the integrity of the syringe and needle has not been compromised and check the expiry date prior to use.
- Do not mix Elevate Kiss, Elevate Ultra and Elevate Voluma with other products or inject the products into the blood stream.
- Do not inject Elevate Kiss, Elevate Ultra and Elevate Voluma in the area selected for treatment prior to a permanent implant site or in connection with other on-going or planned skin treatments, e.g. laser treatment, in the same area.
- Do not use Elevate Kiss, Elevate Ultra and Elevate Voluma in patients with infections or inflammations in the intended injection area.
- Please note that injections of Elevate Kiss, Elevate Ultra and Elevate Voluma can reactivate a latent herpes simplex virus or cause angioedema-type swelling.
- If using tumescent, this may have an effect on the platelet function, such as anti-inflammatory drugs or aspirin, the patient may to a larger extent experience bruising or bleeding at injection sites.
- Elevate Kiss, Elevate Ultra and Elevate Voluma have not been tested on pregnant or breast-feeding women or on children.

CONTRAINDICATIONS

- Elevate Kiss, Elevate Ultra and Elevate Voluma must not be administered to patients with any acute or chronic skin disease or inflammation (such as pimples, rashes or hives) within or close to the area selected for treatment.
- Elevate Kiss, Elevate Ultra and Elevate Voluma must not be administered to patients with bleeding disorders or in patients who are taking thrombolytics or anticoagulants.

PATIENT INFORMATION

- After treatment with Elevate Kiss, Elevate Ultra and Elevate Voluma, undesirable short-term reactions may occur at the injection site, such as pain, redness, swelling, bruising and/or itching. The reactions normally appear shortly after injection and may last up to one week after injections into the skin and up to 2 weeks after lip injections. These types of reactions are common with all types of injections. If the reactions persist or become discomfort or if other reactions occur, such as oedema or prolonged swelling, ante rear the injection site, contact the physician who administered the treatment.
- Elevate Kiss, Elevate Ultra and Elevate Voluma break down naturally over time and are resorbed by the surrounding tissue. Durability varies from patient to patient and the site of the implant.
- The patient must not expose the treated area to extreme temperatures (e.g. sun bathing and solarium) during the first week after treatment or until any undesirable reactions have fully healed. Strenuous exercise should be avoided within the first two days following treatment.
- The patient should be informed not to apply make-up prior to, and up to at least 12 hours after, treatment.

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General directions

- Strict aseptic technique must be performed at all times.
- Elevate Kiss and Elevate Ultra should be injected intradermally under the wrinkle or fold, increasing volume and adding lift and shape.
- When correcting lip lines, Elevate Kiss and Elevate Ultra should be injected in the vermilion border.
- Elevate Voluma should be injected into the subcutaneous tissue or supraperiosteal, increasing volume and adding lift and shape.
- Elevate Kiss, Elevate Ultra and Elevate Voluma are supplied with two sterile needles. The needle size is dependent on the viscosity of the product and the skin thickness.
- Injection is easier to prevent infections is always to be exercised.
- Do not inject Elevate Kiss, Elevate Ultra and Elevate Voluma too deep as this will shorten the effect of the treatment.
- Do not inject Elevate Kiss, Elevate Ultra and Elevate Voluma too superficially as this may result in visible lumps and/or discoloration of the skin.
- Inject slowly into the dermis while pulling the needle backwards.
- Complete the injection a couple of millimetres before the entry point in the skin to prevent material from leaking out of the injection site.
- Correction should only be made to 100% of the desired volume effect. Do not overcorrect.
- A new needle should be used when switching to a new treatment site.
- For each treatment site a maximum dosage of 2 ml is recommended.
- Following gel injection, the product may be assessed for each patient individually. Patients with unavailable expectations should not be subject to treatment.
- The patient's medical background should always be checked before treatment in order to ensure that there are no contraindications. The patient should be informed about the indications, contraindications as well as warnings and precautions.
- The need for any pain relief should be determined.
- In order to achieve and sustain the maximum desired result, the patient should be informed of the need for additional treatment(s).

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- For information on how to assemble the needle and syringe, please see below.
- To avoid the risk of infection during implant, ensure that the skin has been completely disinfected with alcohol, or other antiseptic solution. Standard practice

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
Elevate® Strong

PREPARATION

Contents	Elevate Strong
Sodium hyaluronate, cross-linked	20 mg/ml
Sodium chloride solution	q.s.
Physiological buffer	q.s.

DESCRIPTION

Elevate Strong is a transparent, injectable implant in the form of a viscous, sterile, non-pyrogenic gel. Elevate Strong is composed of cross-linked sodium hyaluronate of non-animal origin.

Elevate Strong is supplied in a graded, disposable, 1.0 ml glass syringe. The product is for single use only. Two disposable sterile needles are provided with each syringe. Information about the needle size and sterilization method can be found on the outer package of the needle. The volume in each syringe is stated on the outer package. Patient record labels are included in each package to ensure product traceability.

INDICATIONS

Elevate Strong is intended for treatment of medium to deep wrinkles and folds, and lip enhancement. The product should be injected intradermally.

WARNINGS AND PRECAUTIONS

- Always observe the general precautions applicable to intradermal injections.
- Elevate Strong contains traces of BDDE (1,4-butanediol diglycidyl) ether, less than 0.00011% used as cross-linking agent, which in rare cases can cause irritation to the skin.
- Elevate Strong may only be used as an intradermal implant by properly trained, qualified physicians or authorized staff members in accordance with national legislation.
- The customary intradermal injection technique must be used during treatment.
- Ensure that the integrity of the syringe and needle has not been compromised and check the expiry date prior to use.
- Do not mix Elevate Strong with other products or inject the product into the blood stream.
- Do not inject Elevate Strong directly adjacent to a permanent implant site or in connection with other on-going or planned skin treatments, e.g. laser treatment, in the same area.
- Do not use Elevate Strong in patients with infections or inflammations in the intended injection area.
- Please note that injections of Elevate Strong can reactivate a latent herpes simplex virus or cause angioedema-type swelling.
- If using substances that may have an effect on the platelet function, such as anti-inflammatory drugs or aspirin, the patient may to a larger extent experience bruising or bleeding at injection sites.
- Elevate Strong has been tested on pregnant or breast-feeding women or on children.

CONTRAINDICATIONS

- Elevate Strong must not be administered to patients with any acute or chronic skin disease or inflammation (such as pimples, rashes or hives) within or close to the area selected for correction.
- Elevate Strong must not be administered to patients with bleeding disorders or in patients who are taking thrombolytics or anticoagulants.

PATIENT INFORMATION

- After treatment with Elevate Strong, undesirable short-term reactions may occur at the injection site, such as pain, tenderness, redness or swelling. The reactions normally appear shortly after treatment and may last up to one week after injections into the skin and up to 2 weeks after lip injections. These types of reactions are common with all types of injections. If the patient experiences severe discomfort or if other conditions, such as oedema or prolonged swelling, arise near the injection site, contact the physician administered the treatment.
- Elevate Strong breaks down naturally over time and are reabsorbed by the surrounding tissue. Durability varies from patient to patient and the site of the implant.
- The patient must not expose the treated area to extreme temperatures (e.g. sun bathing and solarium) during the first week after treatment or until any undesirable reactions have fully healed. Strenuous exercise should be avoided within the first two days following treatment.
- The patient should be informed not to apply make-up prior to, and up to at least 12 hours after, treatment.

DIRECTIONS FOR USE

N.B. DO NOT USE THE PRODUCT IF THE PACKAGE HAS BEEN OPENED OR DAMAGED!

General directions:

- Strict aseptic technique must be performed at all times.
- Elevate Strong should be injected intradermally under the wrinkle or fold, increasing volume and adding lift and shape. When correcting wrinkles and folds, although they are smoothed out, the line will always remain.
- The result of the treatment is based on using a correct injection technique. Elevate Strong must therefore only be administered by well-trained professionals in accordance with local legislation.
- The expected results of treatment must be assessed for each patient individually. Patients with unrealizable expectations should not be subject to treatment.
- The patient's medical background should always be checked before treatment in order to ensure that there are no contraindications. The patient should be informed about the indications, contraindications as well as warnings and precautions.
- The need for any pain relief should be determined.
- In order to achieve and sustain the maximum desired result, the patient should be informed of the need for additional treatment(s).

Directions before and during treatment:

- For information on how to assemble the needle and syringe, please see below.
- To avoid the risk of infection during implant, ensure that the skin has been completely disinfected with alcohol, or other antiseptic solution. Standard practice to prevent infections is always to be exercised.
- Do not inject Elevate Strong too deep as this will shorten the effect of the treatment.
- Do not inject Elevate Strong too superficially as this may result in visible lumps and/or discoloration of the skin.
- Inject slowly into the dermis while pulling the needle backwards.
- Complete the injection a couple of millimetres before the entry point in the skin to prevent material from leaking out of the injection site.
- Correction should only be made to 100% of the desired volume effect. Do not overcorrect.
- A new needle 2ml shall be used when switching to a new treatment site.
- For each treatment site a maximum dosage of 2 ml is recommended.
- Following each injection, massage the treated area gently but thoroughly to ensure the product to follow the contours of the surrounding tissue.
- If stored below room temperature, allow product to equilibrate to room temperature for at least 20 minutes prior to use.

Disposal and single use:

- The syringe, needle and any unused material must be discarded immediately after the treatment.
- No material remaining after treatment may be reused for more than one patient, be saved or be re-sterilized.
- Disposal should be in accordance with accepted medical practice and applicable national or local laws and guidelines.

N.B. ANY ADVERSE EVENTS MUST BE REPORTED TO THE LOCAL DISTRIBUTOR OR DIRECTLY TO THE MANUFACTURER.

ASSEMBLY OF NEEDLE AND SYRINGE

- To avoid separation of the needle and syringe during implant, please note the following instructions:
- Remove the tip cap as displayed (pictures 1 and 2).
- Hold the syringe by using the thumb and forefinger (3) to hold firmly around both the glass syringe barrel and the luer-lock adapter, then grasp the needle (still inside the needle protector) with the other hand (3).
- To facilitate proper assembly, both push and rotate firmly at the same time (4).
- Remove the needle protector in a straight, forward motion (5).
- Strict aseptic technique must be performed at all times.

STORAGE AND SHELF LIFE

- Store the products at a temperature of 2-25°C (36-77°F).
- The products must not be frozen.
- Keep out of direct sunlight.
- The expiry date of the product is indicated on the package.

Manufactured by:

Bohus BioTech AB
Trädgårdsgatan 4
SE-452 31 Strömstad, Sweden
E-mail: info@bohusbiotech.com
Web: www.bohusbiotech.com

Distributed by:

Quttainah Medical Group
KIPCO Tower, Floor 45, Office 2,
Al Shuhada St., Kuwait City
Email: info@qmg-kuwait.com
Web: www.qmg-kuwait.com

Elevate®
STRONG
Cross-linked sodium hyaluronate 20 mg/ml

SYMBOLS

Consult instructions for use	
Single-use only	
Steam sterilized product	STERILE STEAM
Sterile / EO needle	STERILE EO
Store at 2-25° C	2°C 25°C
Keep dry	
Do not use if package is damaged	
Keep away from sunlight	
Batch code	LOT
Use by date	
Manufacturer	

MANUFACTURED BY
Bohus BioTech AB, Trädgårdsgatan 4,
SE-452 31 Strömstad, Sweden
Tel: +46(0)526-146 86 Fax: +46(0)526-126 22
Website: www.bohusbiotech.com
E-mail: info@bohusbiotech.com

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